

## TECHNICAL DATA SHEET

# KI-OSE TUBAX

Code AR0006143

**High Efficacy Disinfectant Cleaner for all surfaces efficient on tubercle bacilli**



### DESCRIPTION

Ki-ose Tubax is a ready to use, all-purpose disinfectant and cleaner for surfaces. It has been specially developed for the disinfection and cleaning of surfaces in aircraft. It is used by cabin attendants for many applications in galleys (surfaces in contact with food), and all surfaces on trolleys, in lavatories and in the cabin. Ki-ose Tubax is an effective disinfectant with a broad spectrum included *Mycobacterium Tuberculosis*, which will not harm or damage fabrics, plastic, metals, rubber or other surfaces.

### FEATURES & BENEFITS

- Non-flammable easy-to-use
- Wide spectrum surface disinfectant
- Suitable for use in food preparation areas
- Very light floral fragrance
- Kills 99.99% of most of harmful bacteria, viruses & fungi (see efficacy tests below)
- Active on *Mycobacterium Tuberculosis*



### DIRECTIONS FOR USE

Spray Ki-ose Tubax on the surface to be disinfected/cleaned. Allow 5 minutes of contact times before wiping with a lint free cloth. Do not rinse except if the surfaces will be in contact with food.

*Use biocides with caution. Before use, read the label and information concerning the product.*

### PHYSICAL PROPERTIES

Active ingredient:	0,2% w/w	Didecyldimethylammonium chloride	CAS 7173-51-5
	0.375% w/w	Glutaraldehyde	CAS 111-30-8
Physical State:	Liquid		
Colour:	Colourless to slightly yellow		
Odor:	Light floral perfume		
Shelf life:	1 year		

*Expiry date and batch number printed on each bottle.*

## TECHNICAL DATA SHEET

### EFFICACY :

<b>BACTERICIDAL</b>	
EN 1276	<i>Pseudomonas aeruginosa, Escherichia coli, Staphylococcus aureus, Enterococcus hirae</i> surrogated bacteria for Enterobacteria
EN 1040	<i>Pseudomonas aeruginosa, Staphylococcus aureus</i>
EN 13623	<i>Legionella pneumophila</i>
EN 13697	<i>Pseudomonas aeruginosa, Escherichia coli, Staphylococcus aureus, Enterococcus hirae</i> surrogated bacteria for Enterobacteria
<b>FUNGICIDAL</b>	
EN 1275	<i>Candida albicans</i> , surrogated fungus for <i>Aspergillus Niger</i>
EN 13697	<i>Aspergillus brasiliensis, Candida albicans</i> , surrogated fungus for <i>Aspergillus Niger</i>
<b>VIRUCIDAL</b>	
EN 14476	Influenza A (H1N1) surrogated virus for lipophilic viruses (Ebola, Coronavirus, Flu, Hepatitis, HIV, Rotavirus), murine norovirus, adenovirus
EN 16777	Adenovirus, murine norovirus
<b>MYCOBACTERICIDAL</b>	
NF T72-151	<i>Mycobacterium Tuberculosis</i>
NF T72-301	<i>Mycobacterium Tuberculosis</i>

### CERTIFICATIONS AND APPROVALS

- Complies with AMS1453 and Boeing D6-7127 standards
- CML Airbus
- Authorized for food contact in accordance with the French decree dated 8th September 1999 concerning the cleaning products for equipment's that may come into contact with foodstuffs.
- Registered with the French Ministry for Ecology, Sustainable Development and Energy (MEDDE) - Inventory N° 64501 / product type 2 and 4 – As part of EU Biocidal Products Regulation N°528/2012.

### ORDERING INFORMATION

Item Code	Packaging	Units/carton	Carton weight	Carton dimensions (cm)
AR0006143	5L drum	4 units	20 kg	38 x 28 x 30

### STORAGE & DISPOSAL

Store in a cool, dry, ventilated place. Keep out of reach of children. Dispose of container in accordance with local regulations.

**WARRANTY** – All statements, information and data presented herein are believed to be accurate and reliable but are not to be taken as a guarantee, expressed or implied, for which seller assumes legal responsibility and they are offered solely for your consideration, investigation and verification. Statements or suggestions concerning possible use of this product are made without representation or warranty that any such use is free of patent infringement and are not recommendations to infringe on any patent. Created 22 April 2020 Modified 27 October 2020 Date Printed 17/12/2020 4:25 PM